


*tie*

*TOURNIQUET  
IN  
EMERGENCY*



Working together

- a) Adjustment of windlass rod every 90°.
- b) Compression plate of 7.4 cm length for use by children and adults.
- c) Resistant windlass fixing clip.
- d) Blocking tab for windlass rod.
- e) Rigid and unbreakable windlass rod.




## Velcro® Self-Engaging Strap

- Maximum grip guaranteed in both horizontal and vertical tractions.
- Grip surface of 4.2 cm width: Up to 20 % wider than other tourniquets.

## Double Safety

- Allowing to allocate remaining band via buckle.
- Preventing any likely slackers in tourniquet use.



## Double Routing Buckle

- Improved in fastening in arms (1) and legs (2).

Weight: 66 g ■ Length: 15 cm (folded) 90 cm (unfolded)

ref. 0200130 BLACK (NSN 6515332153184)

ref. 0200131 NATO GREEN (NSN 6515332077128)

## Tourniquet pouch

- Tourniquet pouch including a hook-n-loop adjustable front closure with red tip indicating opening direction and a single rear MOLLE strap to easily attach onto belts, bags and chest rigs.

ref. 0200132



Puerto de Navacerrada 3  
28935 Móstoles, Madrid, SPAIN  
Tel.: +34 91 616 6000 • Fax: + 34 91 616 4892  
export@hersill.com  
[www.hersill.com](http://www.hersill.com)

